

# Device and Ritual Collection

## Instructions



### Ritual Heat

---

Stir the mask for up to 5 minutes until it reaches an even, thick texture. Then, apply a thin layer to the forehead, nose, chin, cheeks, neck, and upper chest. Leave it on for 10 minutes.

(If purchased without the Treat Serum), activate the mask by mixing it with a small amount of water and letting it rest for 10 minutes.

Once complete, wipe down the face with warm water and pat the skin dry.

### Ritual Treat

---

Apply 2 - 3 pumps on to fingertips & gently begin to massage the serum on to every area the mask is covering. When serum begins to feel warm it will indicate effective application. leave it on for 5 minutes.

### Elevare Red Device

---

You may now begin to use your *Elevare Red Device* for up to 5 minutes per area. Lightly gliding over all surface areas of application.

After completion, remove the excess treatment with a warm towel and clean the silicone cap with warm water and soap. Once fully removed, wash your face with cold water and apply the final step.

### Ritual Complete

---

Apply a generous amount to the face and neck before bed.